




Coned Plug


Features:

- For blocking the flow of liquids and gases
- Standard threads and curvature, no gap, no wobble and good sealing with the corresponding threads
- Imported raw material, non-toxic, tasteless, heat and corrosion resistant.
- PP material is widely used in biomedical, environmental protection and food industries; PTFE material is also suitable for chemical industry and other fields


Unit: mm



WDD-28P



WDD-28E



WDD-06P/E

Model No.	Thread	Material
Coned Plug		
WDD-28P	1/4-28UNF	PP natural
WDD-28E	1/4-28UNF	PTFE natural
WDD-06P	M6	PP natural
WDD-06E	M6	PTFE natural